

SPEC SHEET

Softgel Encapsulation Machine

RG0.8-110, 200, 250, 300C MODELS

INTRODUCTION

A Softgel encapsulation machine fills completely airtight sealed gelatin shells (that exist in one piece) with semi-solids or liquids. Captek softgel encapsulation machines work to the highest specification and reliability levels with unrivaled installation and after sales support.

KEYWORDS Softgel | encapsulation | gelatin | capsules





ATTRIBUTES

- Large touch screen with fault diagnostic functions, operating parameter presets
- Premium aviation aluminum alloy construction
- Fast, stable pneumatic mold pressing
- Independent, variable drive motor, easy to adjust, low noise and smooth operation
- Gelatin sheet cooling mechanism - low energy consumption and high refrigeration efficiency
- Net hopper for soft capsule discharge
- The material hopper can adjust temperature, mixing speed for different product requirements.
- Functional design for easy operation and maintenance

TECHNICAL PARAMETERS

Model	RG0.8-110C	RG2-200C	RG2-250C	RG2-300C
Roller Die Rated Speed	0~7r/min	0~7r/min	0~5r/min	0~5r/min
Roller Die Size	Φ72×110mm	Φ103×200mm	Φ150×250mm	Φ150×300mm
Single Piston Feeding Volume	0~0.8ml	0~2ml	0~2ml	0~2ml
Dimensions	1500×860×1700 mm	1990×1040×2100 mm	2420×1200×2200 mm	2420×1250×2200 mm
Weight	700kg	1000kg	1800kg	1900kg
Total Power	3.13kw	4.75kw	7.55kw	8.15kw

Machine Specs (All Models)	Power Supply	380V; 50Hz
	Noise	<75dBA



While the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.

CONNECT WITH US

+86 577 6688 6278 | info@capsulatech.com | www.captteksolutions.com